Displaying Image Information

John Muschelli

Johns Hopkins Bloomberg School of Public Health

January 25, 2015

Reading In Data


The readNIfTI command (oro.nifti package) can read in NIfTI file (compressed or not) into a nifti object. We will read in the NIfTI (T1-weighted) image from the previous session:

```
library(oro.nifti)
print({nii = readNIfTI(fname = "Output_3D_File")})
NTfTT-1 format
            : nifti
 Type
 Data Type : 4 (INT16)
 Bits per Pixel : 16
 Slice Code : 0 (Unknown)
 Intent Code : 0 (None)
 Qform Code : 2 (Aligned_Anat)
  Sform Code : 2 (Aligned_Anat)
  Dimension : 512 x 512 x 22
 Pixel Dimension: 0.47 x 0.47 x 5
 Voxel Units : mm
 Time Units : sec
```

Visualizing a Slice

The nifti object is a 3D array (see ?array) with header information. We can use the image function (graphics package) to visualize a slice (slice 20 in the z-direction/axial):


image(nii[,,20])



Visualizing a Slice

graphics::image uses heat.colors(12) for coloring, which is not useful for this task. We can either set the colors manually, or use the oro.nifti::image function. The function is still image, but we don't pass in a slice, but the nifti object and specify the slice z=20:


```
image(nii, z = 20, plot.type='single')
```



Visualizing a Slice

If plot.type is not 'single', image.nifti defaults to plotting ALL slices with data, even if z is specified (also called a "lightbox"):

image(nii, z = 20)



Visualizing all 3 planes

To show all 3 planes (axial, sagittal, and coronal) of an image, we can use the orthographic function:

orthographic(nii)



Histograms

What about the **data**? We can do normal operations, such as histograms of the image intensities and intensities over 20:

```
par(mfrow=c(1,2));
hist(nii, breaks = 2000); hist(nii[nii > 20], breaks = 2000
```




Image Overlays

We can do overlays as well, where we have one image and color it by a second. For example, we plot slice 10 and highlight values between 300 and 400 (next slide we discuss the code):





Image Overlays: Explained

We load the fslr package, which has helper functions for nifti objects. The nii > 300 & nii < 400 operation returns an array, not nifti object. The niftiarr command takes in a nifti object and array and returns a nifti object with the array in the data slot. We then set any 0 in mask to NA so the overlay (oro.nifti package) will not mask out data.

Image Overlays: 3 Planes

We can perform the same operation of overlaying, but in all 3 planes:



Functions discussed here

- readNIfTI: read in data
- graphics::image: display matrix data
- oro.nifti::image: display nifti data
- oro.nifti::orthographic: display 3-planes of an image
- oro.nifti::overlay: display overlay of 2 images (NA are not plotting in y image)
- fslr::niftiarr(x,y: Copy nifti header from object x and put in new array y into data part